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EN Percutaneous Introducer Set - Instructions for use

Intended use

The **Introducer sets** of PHS Medical GmbH are sterile accessories used exclusively for the insertion of the intravascular catheters included in the subcutaneously implantable port systems of PHS Medical GmbH.

The product is used to assist the physician in implanting the PHS Medical port systems through a minimally invasive procedure known as the Seldinger technique. The duration of use depends on the duration of the surgery and can range from five to 30 minutes. Different **introducer sets** are available. The appropriate set is selected by the attending physician depending on the catheter size required.

(see table below)

Contents of the percutaneous introducer set

- Puncture cannula (Seldinger cannula)
- 10ml syringe with Luer-Slip
- Guidewire with dispenser and insertion aid
- Split-sheath with dilator
- Tunneler

Storage conditions

The product must be stored clean, dust-protected, dry and free from pests. The product must be stored in such a way that it is protected from damage and direct sunlight. The storage temperature should be between 5 and 35 °C (+41 °F to +95 °F).

Indication

The **Introducer sets** from PHS Medical GmbH are medical devices that can be used, if required, for the implantation of port systems distributed by PHS Medical GmbH. Portsystems provide permanent access to the central venous vascular system and are used, for example, for the application of vaso-stimulant drugs such as chemotherapeutic agents, for long-term infusion therapy or for parenteral nutrition.


Contraindications

The use of **Introducer Sets** may be contraindicated under the following conditions:

- Presence of severe chronic obstructive pulmonary disease
- Presence or suspected presence of vascular occlusion
- Thoracic irradiation
- Known infection, bacteraemia or septicemia
- Presence or suspicion of allergies or intolerance to materials used in the procedure
- Previous venous thrombosis or surgical procedures at or near the site of intended use

WARNINGS:

- This procedure may only be performed by physicians who have been thoroughly trained in it.
- The Tyvek bag containing the tray with the components is a sterile barrier system. If it is damaged or opened, sterility is no longer ensured. The outer folding box is an additional protection for this sterile barrier system.
- For single use on patients. DO NOT REUSE.
- DO NOT resterilize.
- After use, dispose of the product in accordance with your facility's regulations for the disposal of potentially biohazardous materials.
- During venipuncture, care must be taken to prevent air from entering the vein (e.g., Trendelenburg position for air embolism prophylaxis).
- „Pinch-Off“: It has been proven and widely described in the literature that incorrect placement of the catheter in the so-called costoclavicular angle can lead to a „pinch-off“ of the catheter. Therefore, use caution when placing the catheter in the subclavian vein. Make sure that the catheter course outside the vein is not in the area of the intersection of the clavicle and the first rib. (Fig. 15)
- The guidewire must be checked before use. In case of deformation, discard it. DO NOT attempt to straighten a kinked or bent guidewire.
- Never make an attempt to straighten the redresser during the procedure. It must be discarded. Use a new guidewire for the procedure.
- The product must not be modified.
- This product is incompatible with MRI examinations, as overheating or uncontrolled movement of the guidewire may occur.
- Solutions containing alcohol or acetone must NOT be used. They may affect the properties of the plastic parts and cause damage to the product.

Instructions for Use

- In order to check the patency and improve the gliding properties of the guidewire, flush the puncture cannula, split sheath and guidewire with sterile saline solution via the respective luer lock connections using the syringe provided.
- Ideally, place the patient in the head down position with the face facing away from the intended venipuncture site.
- Preparation of the puncture site (cleaning, skin disinfection).
- Cover the area around the puncture site with a perforated surgical drape.
- Fill the enclosed syringe 1/3 to approx. 1/2 full with sterile physiological saline solution, vent and attach the puncture needle.
- Puncture the vein. (Fig. 1)

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Caution! If a vessel is iatrogenically injured or punctured, withdraw the cannula and apply manual pressure to the site for several minutes if necessary. If the pleural space is punctured iatrogenically, withdraw the cannula and check for pneumothorax or hemothorax.

- Aspirate blood to verify the success of the puncture. (Fig. 2)
- Remove the syringe from the puncture needle. (Fig. 3)
- Remove the protective cap from the guidewire insertion aid and discard it. Now pull back the lip protruding from the insertion aid flush to the disc of the insertion aid. (Fig. 4); Attach the insertion aid to the Luer hub of the puncture cannula and then insert the guidewire through the cannula into the vein to the target site. (Fig. 5); Remove the insertion aid and dispenser. (Fig. 6)
- Caution!** If you encounter resistance when advancing or retracting the guidewire, determine the cause by X-ray imaging, before the wire is further advanced or retracted. If the guidewire if significant resistance is felt. If the procedure must be aborted, first remove the puncture needle by pulling it back over the wire before the wire can be retracted.
- Caution!** Since guidewires are sensitive and fragile due to their design, they must be handled with particular care.
- Secure the guidewire and carefully pull the puncture needle back over the guidewire and remove it. (Fig. 7)
- Caution!** NEVER pull the guidewire back over the positioned puncture needle! This may damage the guidewire and enfeather. (Fig. 7)
- In some circumstances, a small incision at the puncture site may be helpful for the subsequent step (skin passage of the split sheath). (Fig. 8)

Caution!

Caution! Check the locking of the dilator with the split-sheath (45° bayonet lock)!

- Slide the split-sheath with the dilator over the guidewire. Then insert it through the skin and into the vessel at an angle of approx. 10° -45° to the skin surface with a slight twisting motion. This twisting motion makes it easier and smoother to insert the split sheath with dilator into the skin. (Fig. 9)
- Attach a clamp to the proximal end of the guidewire prevents the guidewire from being inadvertently fully inserted into the patient.
- Hold the split-sheath, pull back and remove the guidewire, then unlock, pull out and remove the dilator. (Fig. 10)

Caution! Close the opening of the split-sheath with a fingertip to prevent blood loss or air aspiration. (Fig. 11)

- Advance the catheter through the split-sheath into the vein. (Fig. 12)

Caution! When inserting the catheter via the subclavian vein, observe the note „Pinch-Off“ in the „Warnings“ section.

- After placing the catheter, verify the correct position of the catheter tip by fluoroscopy or other suitable means
- Grasp the two handles of the split-sheath with both hands. Now break them away from you to split the lock. (Fig. 13)
- Further split the sheath by pulling the handle wings apart and remove the sheath while pulling it out at the same time. (Fig. 14)

Use of the Tunneler:

- Connect the catheter to the barbed end of the tunneler.
- Tunnel the catheter subcutaneously to the target site.
- Caution!** Avoid damaging skin, bones and surrounding organs!
- Caution!** Avoid kinks, bruises and injuries within along the later used catheter section.
- Disconnect the catheter from the tunneler and follow the instructions for using the port system.

Residual risks and undesirable effects

- Damage to the guidewire
- Premature breakage or kinking of the split-sheath
- Air embolism
- Embolism through the medical device
- Occlusion of the medical device
- Dislocation of the medical device
- Hematoma formation
- Pneumothorax
- Hemothorax and/or brachial plexus injury
- Pulmonary embolism
- Sepsis
- Vascular trauma
- Thrombosis

Note: Please report serious incidents to PHS Medical as well as to the responsible authorities!

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