
PHS MEDiCAL®

C-Port® (CTKP)

C-Port® (SPB)

PICC-NuPort®-CT (LPP)

NuPort®-LP-CT (LPA)

NuPort®-CT (CTP)

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EN Instructions for Use

Intended Use

NuPort®/C-Port® port systems are subcutaneously implantable medical devices used when repeated access to the vascular system is indicated as part of therapy. These products provide the user with an easy-to-find site for cannula insertion to administer fluids, such as drugs, into the vasculature or to draw blood samples. The port is accessed by percutaneous puncture with a (non-coring) port cannula.

The **NuPort®/C-Port®** port systems are available as standard vascular access ports, which are NOT suitable for high pressure injections, and as port systems for CECT high pressure injections. High pressure injections are performed using an appropriate infusion set with a (non-coring) port cannula. Port systems are intended for single use and are explanted and disposed of at the end of therapy.

Caution!

- This product may only be used after comprehensive training by physicians or on their behalf.
- **DO NOT REUSE! DO NOT RESTERILIZE!** Failure to observe these warnings may result in an allergic reaction or infection/cross-contamination in the patient.
- Alcohol, acetone or other solutions containing these substances must **NOT** be used. These solutions can influence properties of the plastic components and lead to their degradation.

Description

The ports of the **NuPort®/C-Port®** family are MRI compatible and consist of an injection septum (silicone) mounted on a titanium or plastic base (PSU) and a silicone or polyurethane radiopaque outlet catheter. This injection septum is designed to withstand multiple punctures with a (non-coring) port cannula while maintaining tightness.

Always use (non-coring) port cannulas supplied with each product to puncture the septum.

For a complete listing of the components included in the port sets and their materials, please refer to the SSPC (see QR code on the back of these instructions for use).

Caution!

DO NOT use standard injection cannulas, as these will cause the septum to leak prematurely.

The septum has a limited service life. Tests have shown that 1,000 punctures/cm² are possible when using a 22-gauge Huber injection cannula. After this number of punctures, the septum may leak. These values may vary individually depending on the size of the cannula used and the technique of drug delivery. In the interest of patient safety, the possibility of leakage and elective replacement should be considered after 1,000 punctures.

Attention!

Always give the patient the enclosed implant ID card and make sure that the patient label with the lot no. and UDI has been applied. This is the only way to ensure traceability. There are also important instructions for use on the ID card.

Please also give the patient the instructions for use, as these should be available to the doctor providing further treatment as well as to the nursing staff.

Accessories

Introducer Sets (Percutaneous Introducers)

(see table „Introducer Set“)

Applicable cannula sizes

The **NuPort®-CT**, **C-Port®**, **C-Port®-CT** and **C-Port®-II** products can be used with (non-coring) port cannulas in sizes **19 G**, **20 G** and **22 G**.

Caution!

With the **PICC-NuPort®-CT** and **NuPort®-LP-CT**, **ONLY** (non-coring) port cannulas in sizes **20 G** and **22 G** may be used for standard injection and **ONLY** (non-coring) port cannulas in size **22 G** may be used for high-pressure injection.

Storage Conditions

The product must be stored clean, dust-protected, dry and free from pests.

The product must be stored in such a way that it is protected from damage and direct sunlight. The storage temperature should be between 5 °C and 35 °C (+41 °F to +95 °F).

MR Compatibility

Conditionally MR safe

The ports of the **NuPort®/C-Port®** family have been tested for safety in an MRI environment at the University of Southern California Hospital MRI Center in accordance with the American Society for Testing and Materials (ASTM). The ports of the **NuPort®/C-Port®** family were found to be **conditionally MR safe**⁶¹⁷⁾.

The products can thus be safely subjected to an MRI scan immediately after placement under the following conditions:

- Static magnetic field of 3 Tesla or less
- Maximum spatial gradient magnetic field of 12,000 Gauss/cm (extrapolated) or less
- Maximum whole-body averaged specific absorption rate (SAR) of 4 W/kg for 15-minute scanning (i.e., per pulse train) specified for the MR system
- First level controlled operating mode for operation of the MRI system.

Attention!

The MR compatibility refers only to the port with catheter!

MR-Induced Heating

In non-clinical testing, the ports* of the **NuPort®/C-Port®** family produced the following maximum temperature rise during MRI examinations performed with a 15-minute scan (i.e., per pulse train) in a 3 Tesla MR system (3 Tesla/128 MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI):

- Whole-body averaged SAR reported for the MR system: 2,9 W/kg
- Calorimetry readings, whole-body averaged SAR: 2,7 W/kg
- Highest temperature change: ≤ 2,4 °C*
- Temperature scaled to a whole-body averaged SAR of 4 W/kg: ≤ 3,6 °C*
- Local SAR at the reference probe position: 3,2 W/kg

(*Based on worst-case C-port DL)

MR Artifacts

The quality of the MR image may be affected if the relevant area is at the position of the product or relatively close to it. Therefore, it may be necessary to optimize MR imaging parameters in the presence of the implant.

Indication

The ports of the **NuPort®/C-Port®** family are intended for use in patients who require: repeated access to the vasculature for injections, infusion of drugs, administration of blood or blood products, and/or collection of blood as part of the therapeutic regimen. Port and catheter size should be selected according to the patient's age and stature. The **NuPort®/C-Port®** port systems can be implanted in different areas of the body.

Areas of application include:

- parenteral nutrition
- chemotherapy
- pain therapy
- endocrinological therapy
- substitution of coagulation factors
- contrast medium infusion in diagnostic radiology
- therapy of severe asthma or cardiac arrhythmia
- blood collection and transfusion

Contraindications

The ports of the **NuPort®/C-Port®** family are contraindicated in the following conditions:

- In patients with known or suspected allergic reactions to the materials used to manufacture the product
- In patients with known infection, bacteremia or sepsis
- In patients with body tissue that cannot adequately support the port or catheter
- In patients with chronic obstructive pulmonary disease (COPD)
- In patients with previous venous thrombosis or surgical procedures at or near the site of intended use
- When a potential access site has been irradiated
- When the patient's tissue factors prevent proper stabilization of the device and/or access

Potential Complications

The use of the ports of the **NuPort®/C-Port®** family is associated with risks known to be associated with the insertion or use of an implanted device or indwelling catheter. These risks include, but are not limited to: infection, erosion of skin, hematoma, necrosis, displacement of device, clot formation, catheter fragmentation, embolization, catheter occlusion, cardiac arrhythmias, cardiac puncture, pericardial tamponade, extravasation of drugs, fibrin sheath, endocarditis, implant rejection, port/catheter migration, perforation or rupture of vessels, thoracic duct injury, thromboembolism, thrombosis, superior vena cava syndrome, pneumothorax, spontaneous misplacement or retraction of catheter tip, vascular erosion, catheter or port related sepsis, air embolism, hemorrhage, catheter embolism.

Attention!

PHS Medical GmbH cannot guarantee the compatibility of the drugs used with the port system.

Catheter Fracture (pinch-off)

Pinch-off can occur when the catheter is compressed between the clavicle and the first rib (costoclavicular gap). This complication occurs in approximately 1-5% of implantations.

Clinical signs:

- Difficulty with blood sampling, excessive pressure during infusion of fluids.
 - Patient must be repositioned when infusing fluids or for blood sampling
- If these complications occur, radiologic evaluation should be performed to determine the cause.

Radiologic Assessment:

Evidence from the chest x-ray is used to grade the degree of deformity using the table below:

Grad	Deformity	Recommended treatment
Grad 0	No deformity	No treatment regarding pinch-off; clarify other causes.
Grad I	Deformation present, but lumen intact	Monitor with chest x-rays in 1-3 months to check for pinch-off progression
Grad II	Deformation present with narrowed lumen	Consider removal of catheter
Grad III	Catheter damaged or fractured	Remove catheter immediately

Prevention of pinch-off

- By placement via the internal jugular vein, the risk of pinch-off is almost completely eliminated.
- When placing the catheter in the subclavian and superior vena cava, ensure that the catheter does not pass through the intersection of the clavicle and the first rib¹⁾⁵⁾.
- Catheters inserted via the subclavian vein should be inserted into the vessel at the junction of the lateral and middle thirds of the clavicle. (Do not insert the catheter too far medially).
- X-rays should be used to ensure that the catheter is not trapped between the first rib and the clavicle.

Note: Please report serious incidents to PHS Medical GmbH as well as to the responsible authorities!

Information on Sterility

This product is supplied sterile and non-pyrogenic.

The Tyvek pouch containing the tray with the components is a sterile barrier system. If it is damaged or opened, sterility is no longer assured. The outer folding box is an additional protection of this sterile barrier system.

This product is intended for single use only. Products from PHS Medical GmbH must not be reimplanted under any circumstances. Products that have been contaminated with blood must not be reused or resterilized.

Implantation Instructions

Use the surgical procedure and sterile technique for implantation of the port system that, in your experience, best suits the patient's needs. The following recommendations and safety-related instructions for surgical insertion are intended to assist in the safe and long-term use of the port systems of the **NuPort®/C-Port®** family.

Caution!

The product may **ONLY** be used if the sterile barrier system (Tyvek pouch) is intact.

- 1) Visually inspect the port and catheter before implantation. Do not use if holes, cracks or surface contamination can be found.
- 2) Purge all air from the port and catheter before inserting the product into the patient.
- 3) Place the port in a suitable location. Ensure that there is an adequate amount of adipose tissue between the port and the skin so that the port chamber can be easily palpated. In obese patients, a more parasternal placement of the port chamber may be advantageous. If possible, a site should be chosen that is comfortable and convenient for the patient (e.g., for female patients, the bra strap area should be avoided, etc.). This provides the patient with an inconspicuously implanted vascular access that is not disturbed by normal body movement.

Caution!

The patient should be in the head-down position („Trendelenburg position“) with the face turned away from the intended site of venipuncture to prevent air embolism!

- 4) Insert the catheter into the vein using either a surgical technique („cut-down“) or a percutaneous introducer of appropriate size (see table „Introducer Set“).
- 5) Place the tip of the catheter in an area of high blood flow. It is recommended to check the correct placement of the catheter tip in the superior vena cava (ideally between the 5th and 6th rib) under fluoroscopy.

Attention!

Especially with the **PICC-NuPort®-CT** and **NuPort®-LP-CT** as an arm port, the catheter should **ONLY** be inserted in the upper arm! The port can be implanted in the forearm and the catheter tunneled from the upper arm⁶⁾.

- 6) Ensure that the catheter tip is not serrated or occluded during catheter placement. After the catheter is positioned in the vessel, leave sufficient slack in the catheter between the insertion point and the port. This allows the patient to move without stressing the catheter insertion point or pulling the catheter out of its intended location.
- 7) Position the port pocket so that the suture line is not directly over the port. The suture line can be lateral, medial, and above or below the port body, and the pocket should be large enough to avoid stretching or trauma at the suture line. Do not place the port so deep that it is difficult to find, and do not place it so shallow that erosion or pressure necrosis can be caused. The recommended optimal depth for placement is approximately 5 mm below the skin surface.

Instructions for Attaching the Catheter

The port outlet consists of two components, a special „bulb-shaped“ outlet („olive“) and a catheter lock (Fig. 1). Together, these two components serve to connect the catheter to the port without leakage. They are connected as follows:

- a) Cut the catheter to the appropriate length and wet all components with saline solution for lubrication and easier assembly.
- b) Slide the catheter lock over the catheter (see Fig. 2).
- c) Slide the catheter over the „bulb“ outlet („olive“) so that it is between the last bulb and the port body (Fig. 2).
- d) Push the catheter lock together with the catheter forward until the catheter rests against the port body and the catheter and outlet are completely covered (Fig. 3).
- e) Check the catheter connection by gently „pulling“ on the catheter. (Fig. 4)

Caution!

If the catheter and catheter lock are already connected and then disconnected, the catheter end must be shortened to ensure a safe reconnection.

- 8) Attach the port to the underlying fascia with a minimum of three (3) independent sutures placed around the base. Pass each suture through the holes in the port base and fascia. This will prevent migration of the port and eliminate the risk of the port flipping over in the pocket.

Caution!

During implantation, care must be taken not to constrict the catheter during suturing or to damage it with a needle or scalpel, as damage can cause the catheter to leak or tear quickly.

- 9) After suturing the port to the fascia, flush the pocket with an appropriate antibiotic to

ensure sterility.

- 10) Check patency and flow through the port before wound closure. To do this, first aspirate blood and then flush the port with at least 2 x 10 ml saline solution.

Attention!

Never use syringes smaller than 10 ml syringes.

- 11) Confirm the position of the catheter using X-ray, fluoroscopy, or an imaging modality of your choice.
- 12) Block the port catheter after implantation and each use with a solution recommended by your institution.

Caution!

Some patients may be hypersensitive to heparin or have heparin-induced thrombocytopenia (HIT). In these patients, the use of heparin is contraindicated.

- 13) Dispose of biohazardous materials according to your facility's regulations.

Guidance for Use

Use only sterile pyrogen-free products when injecting, infusing or sampling through the ports of the **NuPort®/C-Port®** family.

Caution!

Your hospital's protocols for care and use of the port may differ from these instructions for use. In this case, the responsible physician at your institution must decide which procedure to use and what constitutes appropriate patient care.

Recommended accessories:

- Sterile gloves and sterile mask
- Accessories for surgical cleaning and care of the injection site
- Sterile (non-coring) port cannula (select size appropriate for intended use)
- Sterile syringe and medication
- Sterile syringe with appropriate solution to block the catheter
- Suitable anesthetic (optional)
- Sterile unidirectional valve or three-way stopcock (optional)

Procedure - Standard Vascular Access

Caution!

- All steps in standard vascular access must be performed under aseptic conditions.
- **DO NOT** use syringes smaller than 10 ml syringes under any circumstances.
- **DO NOT** attempt to perform high pressure injections through a port that is not labeled „CT Compatible / Power Injectable“.
- Observe the note on the cannula sizes to be used under the section „Accessories“!

- 1) Clean the injection site under aseptic conditions using an appropriate surgical cleaning method.
- 2) Using sterile gloves, palpate the implantation site to locate both the port and the septum (Fig. 5).
- 3) To gain access to the port, isolate its top (Fig. 6). Then stabilize the port with your thumb and index finger (Fig. 7) and puncture the skin and septum with the (non-coring) port cannula or a (non-coring) port cannula infusion set. Fill the port cannula or infusion set with saline solution beforehand. Insert the port cannula vertically into the port until it reaches the port base (Fig. 8).

Caution!

During these procedures, the cannula tip must be held securely at the port bottom, under the septum, to avoid injecting the drug into the subcutaneous tissue. Avoid tilting or rotating the cannula and syringe once they are in the septum.

However, the pressure exerted with the cannula tip on the port floor must not be too great, as this will cause bending („fishhooking“) of the cannula tip. Especially with the **PICC-NuPort®-CT** and **NuPort®-LP-CT**, this will lead to excessive damage to the septum, which may cause extravasation. The size of the „fishhook“ depends on the pressure applied by the user during puncture, but is also strongly influenced by the port cannula used (size and design). To avoid excessive damage to the port by puncturing with (non-coring) port cannulas, we recommend:

- avoid excessive pressure on the cannula during access.
 - not to puncture the port at the same point of the septum each time
 - always use a new port cannula when puncturing again (also for puncture correction)
 - to always insert the (non-coring) port cannulas vertically to prevent coring.
- 4) Before each use, perform a cannula position check by slowly aspirating blood. Discard the blood and flush the port with at least 2 x 10 ml of physiological saline solution.

Caution!

- **DO NOT** attempt to inject through a port that shows signs of blockage or pinch-off, as this may result in damage to the port system or the patient.
 - When administering drugs with an infusion pump, the pressure must not exceed 8 bar (116 psi)!
- 5) Perform injections with continuous, firm pressure using a syringe or a suitable infusion pump.
 - 6) Flush the port with at least 10 ml of physiological saline solution after each injection of a drug. Flush the port with at least 2 x 10 ml saline solution after each blood collection or injection of highly viscous fluids, such as parenteral nutrition solution.

Caution!

Flushing is crucial for the continued use of the port, as this reduces the occurrence of catheter occlusions. PHS Medical GmbH recommends the „push-pause technique“ during the entire flushing process.

- 7) If you wish to apply another „catheter block“ than the saline solution you just administered, you may now do so over the port cannula that is still in place, according to your institution's recommendations.

Caution!

Some patients may be hypersensitive to heparin or suffer from heparin-induced thrombocytopenia (HIT). In these patients, the use of heparin is contraindicated.

- 8) Withdraw the port cannula while applying pressure to the plunger of the syringe.

Dispose of used components according to your facility's regulations.

- 9) Examine the injection site closely. Palpate the area to ensure that the injection was delivered to the desired site.

Measures for Mitigating Complications

- Adherence to care standards during port puncture/port cannula exchange
- Application of aseptic technique
- Checking the puncture site for signs of infection
- Questioning the patient about unusual sensation during flushing/infusion
- For the **PICC-NuPort®-CT** and **NuPort®-LP-CT**, observe the information on the cannula sizes to be used under the section „Accessories“!
- Application of a suitable flushing technique
 - o Sufficient flushing volumes
 - o Push-pause technique (pulsed flushing, approx. ½-1-second interval)
 - o Positive pressure during tube clamping / during port cannula removal

Procedure - CECT High Pressure Injection

Only the **NuPort®-CT/C-Port®-CT** port systems labeled „CT Compatible / Power Injectable“ may be used for the administration of contrast medium using a high-pressure injector. Their special design allows high-pressure injection of contrast medium with a high flow rate when a suitable (non-coring) port cannula infusion set is used and the recommended maximum infusion rates and pressure settings are adhered to.

Caution!

- **DO NOT** attempt to perform a high pressure injection through a port that is not labeled „CT Compatible / Power Injectable“ or that shows signs of „pinch-off“ as this may result in damage to the port system or the patient.
- Use only appropriate (non-coring) port cannula infusion sets for high-pressure injection.
- Before performing high-pressure injection, warm the contrast medium to body temperature.
- Before injection, make sure that the contrast medium viscosity is less than or equal to 11.8 mPas.

- 1) Prepare the port for high pressure injection according to steps 1) to 5) described under „Procedure - Standard Vascular Access“.
- 2) Perform the high-pressure injection according to the manufacturer's instructions for use for the high-pressure injector and/or according to your institution's guidelines.
- 3) Upon completion of the procedure, follow steps 7) through 9) under „Procedure - Standard Vascular Access“.

Caution!

When performing high pressure injection through the **NuPort®-CT/C-Port®-CT** port systems, do not exceed the pressure limit of 300 psi (~ 21 bar) or the maximum flow rate setting on the high pressure injector (see table below). Use **ONLY 22 G** cannulas for high pressure injection through the **PICC-NuPort®-CT** and **NuPort®-LP-CT**.

Maximum Flow rates and Pressure*:

	NuPort®-CT/C-Port®-CT			PICC-NuPort®-CT NuPort®-LP-CT
Cannula Size	19 G	20 G	22 G	22 G
Catheter Size	6,6F / 7,5F 8,0F / 9,0F 9,6F	6,6F / 7,5F 8,0F / 9,0F 9,6F	5,0F / 6,6F 7,5F / 8,0F 9,0F / 9,6F	5,0F / 6,6F / 7,5F 8,0F / 9,0F / 9,6F
Max. Flow Rate	5 ml/s	5 ml/s	2 ml/s	2 ml/s
Max. Pressure	300 psi / ~21 bar			300 psi / ~21 bar

*The maximum flow rate may deviate from the specified values depending on the cannula manufacturer! Please observe the corresponding instructions for use!

Blood Collection

When taking a blood sample, the following procedure and associated precautions must be followed:

- 1) Prepare the port for blood collection according to steps 1) to 3) above under „Procedure - Standard Vascular Access“.
- 2) Aspirate 5 ml of blood and discard it with the syringe. If aspiration is not possible, recheck the position of the cannula. Connect syringes as needed to obtain the necessary blood samples.
- 3) After completing the procedure, follow steps 6) to 9) under „Procedure - Standard Vascular Access“.

Attention!

It is extremely important to flush the port with at least 2 x 10 ml saline solution after blood collection. If blood remains in the port system for a longer period of time, occlusion of the system may occur.

Care and Maintenance

PHS Medical GmbH recommends flushing of the ports of the **NuPort®/C-Port®** family with saline solution⁹⁾ every three months during extended periods of inactivity to prevent occlusion. Your hospital's protocols for port care and maintenance may differ from this recommendation. In this case, the physician in charge at your facility must decide how often to flush the port system.

To flush the port, follow steps 1) to 4) and 7) to 9) described under „Procedure - Standard Vascular Access“.

Cleaning Blocked Catheters

If you detect resistance during an injection or infusion, stop the procedure immediately. Excessive pressure can cause damage to the port system or the patient. Resistance can be caused by an occluded catheter or port chamber, as well as catheter entrapment. If catheter entrapment can be ruled out by x-ray, consider clearing with an appropriate lysing agent such as urokinase²⁾³⁾. An initial dose of 5,000 units is recommended, but consult the drug package for instructions and possible complications⁴⁾ when using urokinase.

Explantation

If the port catheter is no longer required for further therapy, it should be explanted in accordance with your institution's regulations and disposed of.

If you encounter resistance when removing the catheter, this may be due to the formation of a fibrin layer around the catheter. In some cases, this can calcify and make it difficult to remove the catheter. Take appropriate measures to remove the catheter without tearing it off.

After explantation, always ensure that the port, catheter fixation, and catheter (using the length markings on the catheter) have been completely removed.

Note: a calcified fibrin sheath can lead to misinterpretation of the X-ray images, as it may look like a piece of catheter remaining in the body. It is therefore essential to ensure that the catheter has been completely removed. If a new port catheter is to be implanted, it must be ensured that the fibrin sheath has been destroyed beforehand so that the new catheter is not inserted into this fibrin sheath.

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PHS Medical Introducer Set	Part No.	Split-Sheath Internal Diameter	Dilator Length	Dilator Internal Diameter	For use with PHS Medical Portsystem
Percutaneous Introducer (6F)	PI-0006	0.083 Inch (2,1 mm)	17,75 cm	0.039 Inch (0,99 mm)	NuPort®/C-Port® with 5F-Catheter
Percutaneous Introducer (7F)	PI-0007	0.096 Inch (2,4 mm)	17,75 cm	0.039 Inch (0,99 mm)	NuPort®/C-Port® with 6F-/6,6F-Catheter
Percutaneous Introducer (8F)	PI-0008	0.109 Inch (2,7 mm)	17,75 cm	0.039 Inch (0,99 mm)	NuPort®/C-Port® with 7,5F-Catheter
Percutaneous Introducer (9F)	PI-0009	0.122 Inch (3,0 mm)	17,75 cm	0.039 Inch (0,99 mm)	NuPort®/C-Port® with 8F-Catheter
Percutaneous Introducer (10F)	PI-0010	0.135 Inch (3,4 mm)	17,75 cm	0.039 Inch (0,99 mm)	NuPort®/C-Port® with 9,0-/9,6F-Catheter
Percutaneous Introducer – Micropuncture Set (5,5F)	PI-055M	0.076 Inch (1,9 mm)	9,25 cm	0.019 Inch (0,48 mm)	NuPort®/C-Port® with 5F-Catheter
Percutaneous Introducer – Micropuncture Set (7F)	PI-007M	0.096 Inch (2,4 mm)	9,25 cm	0.019 Inch (0,48 mm)	NuPort®/C-Port® with 6,6F-Catheter

Catheter / Port Specifications:

Catheter Specifications / Katheterspezifikationen / Especificaciones del catéter Caractéristiques du cathéters / Specifiche del catetere / Especificações do Cateter Характеристики катетера / Προδιαγραφές καθετήρα

Catheter Size	Material	ID (mm)	Priming Volume (mL/10cm)	Gravity Flow Rate* (mL/min)
5.0 F	Polyurethane	1.05	.088	16
6.6 F	Polyurethane	1.28	.113	20
7.5 F	Polyurethane	1.15	.106	18
8.0 F	Silicone	1.40	.154	48
9.0 F	Polyurethane	1.60	.201	52
9.6 F	Silicone	1.58	.196	52

* Flow rate through working length of 20cm / Flußrate durch Arbeitslänge von 20cm / Tasa de flujo a través de longitud de trabajo de 20 cm / Débit à travers une longueur de travail 20 cm / Portata attraverso la lunghezza operativa di 20cm / Taxa de fluxo através do comprimento útil de 20cm/ Скорость потока при учете рабочей длины катетера 20см / Ταχύτητα ροής διαμέσου του ωφέλιμου μήκους των 20 εκατοστών

Port Specifications / Portspezifikationen / Especificaciones del puerto Caractéristiques de la chambre implantable Specifiche della porta vascolare / Especificações do port Характеристики Порта / προδιαγραφές port

Port	Material	Priming Volume	Weight
NuPort-CT	Titanium	0.5 mL	12.2 g
PICC-NuPort-CT	Titanium	0.3 mL	8.0 g
NuPort-LP-CT (LPA)	Titanium	0.3 mL	8.0 g
C-Port-CT	Plastic	0.5 mL	6.3 g
C-Port-II	Plastic	0.5 mL	6.3 g
C-Port	Titanium	0.5 mL	12.5 g

Symbol Erklärung / Symbol Explanation

Symbol	(DE)	(EN)	(ES)	(FR)	(IT)	(PT)	(GR)	(RU)
	Bedingt MR-sicher bis maximal 3-Tesla	Magnetic Resonance (MR) conditional to max. 3-Tesla	Compatible condicional con la RMN hasta un máximo de 3 teslas	Compatibilité avec l'IRM sous certaines conditions, jusqu'à 3 teslas max.	A compatibilità MR condizionata fino a max. 3 tesla	Condicionalmente seguro para RM até um máximo de 3 Tesla	Ασφαλές υπό όρους για μαγνητική τομογραφία έως το μέγιστο 3 Tesla	Условно безопасно при МРТ до 3 Тесла
	Fähigkeit zur Durchführung von Kontrast-mittelinjektionen bis maximal 300 psi bei der CT-Bildgebung	Capable of performing contrast medium injections up to a maximum of 300 psi during CT imaging.	Capacidad para realizar inyecciones de medio de contraste hasta un máximo de 300 psi durante la adquisición de imágenes por TC	Capacité à réaliser des injections de produits de contraste jusqu'à max. 300 psi pour l'imagerie par TDM	Idoneità all'esecuzione di iniezioni di mezzo di contrasto fino a max. 300 psi per imaging TC	Capacidade para efetuar injeções de contraste até um máximo de 300 psi em imagens de TC	Ικανότητα διεξαγωγής ενγύσεων σκιαγραφικού μέσου έως το μέγιστο 21 bar (300 psi) σε απεικόνιση αξονικής τομογραφίας (CT)	Возможность проведения инъекции контрастных веществ под давлением до 300 psi во время КТ

- Alken DR, Minton JP. The "Pinch-Off" Sign: A Warning of Impending Problems with Subclavian Catheters. American Journal of Surgery. 148: 633-636, 1984
- Hurtubise MR, Bottino JC, Lawson M, McCredie KB. Restoring Patency of Occluded Central Venous Catheters. Archives of Surgery, Vol. 115, Feb. 1980.
- Moore CL, Erickson KA, Yanes LB, Franklin M, Gonsalves L. Nursing Care and Management of Venous Access Ports. Oncology Nursing Forum, Vol. 13, No. 3, May/June 1986.
- Fiscus JA, Brown D. Withdrawal Occlusion in VAD. Oncology Nursing Forum, Vol. 16, No. 6, 1989.
- Hinke D.H., Zandi-Stastny D.A., Goodman L.R., et al., "Pinch-Off Syndrome: A Complication of Implantable Subclavian Venous Access Devices." Radiology 177: 353-356, 1990.
- Shellock Frank PhD. Evaluation of Magnetic Field Interactions, Heating, and Artifacts for the C-Port-DL Vascular Access Port. MR Center, Univ. of California, Los Angeles, April 2014.
- Shellock Frank PhD. Evaluation of Magnetic Field Interactions, Heating, and Artifacts for the NuPort-LPA Vascular Access Port. MR Center, Univ. of California, Los Angeles, April 2014.
- JP Götz et al., "Totally Implantable Venous Power Ports for the Forearm and the Chest - Initial Clinical Experience", British Journal of Radiology, 85 (2012) e966-e972
- GA Goossens et al, "Comparing normal saline versus diluted heparin

Hersteller / Manufacturer / Fabricante / Fabricant /
Fabbicante / Fabricante / Κατασκευαστής / Производитель

PHS Medical GmbH
Ederweg 3 · D-34277 Fuldabrück
Germany

Telefon +49 (0)561-998597-0
Email info.phs@numantecgroup.com
www.phs-medical.com

QR SSCP Portsystem

www.phs-medical.de/K07V12F03_P

CE 0483

P/N 01621365 Rev. H 09/25